



OR4C15 Antibody, Biotin conjugated

Product Code	CSB-PA836757OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGM1
Immunogen	Peptide sequence from Human Olfactory receptor 4C15 protein (258-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 4C15 (Olfactory receptor OR11-127) (Olfactory receptor OR11-134), OR4C15
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR4C15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.