



OR4D10 Antibody, FITC conjugated

Product Code	CSB-PA839891OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGI6
Immunogen	Peptide sequence from Human Olfactory receptor 4D10 protein (86-98AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 4D10 (Olfactory receptor OR11-251), OR4D10, OR4D10P
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR4D10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.