




OR5K2 Antibody

Product Code	CSB-PA844083OA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHB8
Immunogen	Peptide sequence from Human Olfactory receptor 5K2 protein (6-24AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 5K2 (Olfactory receptor OR3-9), OR5K2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR5K2

Image



Western Blot
Positive WB detected in: U251 whole cell lysate, U87 whole cell lysate
All lanes: OR5K2 antibody at 1:2000
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 36 kDa
Observed band size: 36 kDa