



OR5T2 Antibody, HRP conjugated

Product Code	CSB-PA850889NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGG2
Immunogen	Peptide sequence from Human Olfactory receptor 5T2 protein (298-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 5T2 (Olfactory receptor OR11-177), OR5T2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR5T2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.