




OR5V1 Antibody

Product Code	CSB-PA871401OA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGF6
Immunogen	Peptide sequence from Human Olfactory receptor 5V1 protein (302-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 5V1 (Hs6M1-21) (Olfactory receptor OR6-26), OR5V1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR5V1

Image



Western Blot

Positive WB detected in: A549 whole cell lysate

All lanes: OR5V1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 37 kDa

Observed band size: 37 kDa