



OR6B3 Antibody, FITC conjugated

Product Code	CSB-PA847734OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGW1
Immunogen	Peptide sequence from Human Olfactory receptor 6B3 protein (259-272AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 6B3 (Olfactory receptor OR2-2), OR6B3, OR6B3P
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR6B3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.