



OR6C1 Antibody, Biotin conjugated

Product Code	CSB-PA850422OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD1
Immunogen	Peptide sequence from Human Olfactory receptor 6C1 protein (259-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 6C1 (OST267), OR6C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR6C1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.