




OSCAR Antibody

Product Code	CSB-PA811628OA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYS5
Immunogen	Peptide sequence from Human Osteoclast-associated immunoglobulin-like receptor protein (265-279AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Osteoclast-associated immunoglobulin-like receptor (Osteoclast-associated receptor) (hOSCAR) (Polymeric immunoglobulin receptor 3) (PIgR-3) (PIgR3) (Poly-Ig receptor 3), OSCAR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OSCAR

Image



Western Blot

Positive WB detected in: THP-1 whole cell lysate, HL60 whole cell lysate

All lanes: OSCAR antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31, 29, 30, 28 kDa

Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.