



PCDHA10 Antibody, Biotin conjugated

Product Code	CSB-PA897505LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5I2
Immunogen	Recombinant Human Protocadherin alpha-10 protein (276-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin alpha-10 (PCDH-alpha-10), PCDHA10, CNRS8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCDHA10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.