



PCDHGB7 Antibody, FITC conjugated

Product Code	CSB-PA896527LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5F8
Immunogen	Recombinant Human Protocadherin gamma-B7 protein (719-929AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin gamma-B7 (PCDH-gamma-B7), PCDHGB7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCDHGB7