



PCDHGB7 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA896527LB11HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y5F8 |
| Immunogen | Recombinant Human Protocadherin gamma-B7 protein (719-929AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Antigen Affinity purified & Affinity purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protocadherin gamma-B7 (PCDH-gamma-B7), PCDHGB7 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | PCDHGB7 |