



PCDHGB7 Antibody, Biotin conjugated

Product Code	CSB-PA896527LD11HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5F8
Immunogen	Recombinant Human Protocadherin gamma-B7 protein (719-929AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin gamma-B7 (PCDH-gamma-B7), PCDHGB7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCDHGB7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.