



PIBF1 Antibody, Biotin conjugated

Product Code	CSB-PA845171LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXW3
Immunogen	Recombinant Human Progesterone-induced-blocking factor 1 protein (1-757AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Progesterone-induced-blocking factor 1 (PIBF) (Centrosomal protein of 90 kDa) (CEP90), PIBF1, C13orf24 PIBF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIBF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.