



POLR2A Antibody, Biotin conjugated

Product Code	CSB-PA02519D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24928
Immunogen	Recombinant Human DNA-directed RNA polymerase II subunit RPB1 protein (76-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerase II subunit RPB1 (RNA polymerase II subunit B1) (EC 2.7.7.6) (DNA-directed RNA polymerase II subunit A) (DNA-directed RNA polymerase III largest subunit) (RNA-directed RNA polymerase II subunit RPB1) (EC 2.7.7.48), POLR2A, POLR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POLR2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.