



PRG4 Antibody, Biotin conjugated

Product Code	CSB-PA018672LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92954
Immunogen	Recombinant Human Proteoglycan 4 protein (25-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteoglycan 4 (Lubricin) (Megakaryocyte-stimulating factor) (Superficial zone proteoglycan) [Cleaved into: Proteoglycan 4 C-terminal part], PRG4, MSF SZP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.