



Ren1 Antibody, Biotin conjugated

Product Code	CSB-PA019561LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08424
Immunogen	Recombinant Rat Ren1 protein (65-264AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Renin (EC 3.4.23.15) (Angiotensinogenase), Ren1, Ren
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Ren1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.