




RPTN Antibody

Product Code	CSB-PA0205010A01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XPR3
Immunogen	Peptide sequence from Human Repetin protein (86-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Repetin, RPTN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPTN

Image



Western Blot

Positive WB detected in: HL60 whole cell lysate, K562 whole cell lysate, HeLa whole cell lysate
All lanes: RPTN antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 91 kDa
Observed band size: 91 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.