



RPTN Antibody, FITC conjugated

Product Code	CSB-PA020501OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XPR3
Immunogen	Peptide sequence from Human Repetin protein (86-103AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Repetin, RPTN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.