



SAMHD1 Antibody, Biotin conjugated

Product Code	CSB-PA020691LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3Z3
Immunogen	Recombinant Human Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 protein (1-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 (dNTPase) (EC 3.1.5.-) (Dendritic cell-derived IFNG-induced protein) (DCIP) (Monocyte protein 5) (MOP-5) (SAM domain and HD domain-containing protein 1), SAMHD1, MOP5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SAMHD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.