



SDHB Antibody, Biotin conjugated

Product Code	CSB-PA020907LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21912
Immunogen	Recombinant Human Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial protein (29-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial (EC 1.3.5.1) (Iron-sulfur subunit of complex II) (lp), SDHB, SDH SDH1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SDHB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.