



SLC1A3 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA021434LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P43003 |
| Immunogen | Recombinant Human Excitatory amino acid transporter 1 protein (146-236AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Excitatory amino acid transporter 1 (Sodium-dependent glutamate/aspartate transporter 1) (GLAST-1) (Solute carrier family 1 member 3), SLC1A3, EAAT1 GLAST GLAST1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | SLC1A3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |