



SLC1A3 Antibody, Biotin conjugated

Product Code	CSB-PA021434LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43003
Immunogen	Recombinant Human Excitatory amino acid transporter 1 protein (146-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Excitatory amino acid transporter 1 (Sodium-dependent glutamate/aspartate transporter 1) (GLAST-1) (Solute carrier family 1 member 3), SLC1A3, EAAT1 GLAST GLAST1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC1A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.