



# SLC1A3 Antibody, Biotin conjugated

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA021434LD01HU  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | P43003  |
| Immunogen           | Recombinant Human Excitatory amino acid transporter 1 protein (146-236AA)   |
| Raised In           | Rabbit  |
| Species Reactivity  | Human   |
| Tested Applications | ELISA   |
| Form                | Liquid  |
| Conjugate           | Biotin  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| Purification Method | >95%, Protein G purified  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Excitatory amino acid transporter 1 (Sodium-dependent glutamate/aspartate transporter 1) (GLAST-1) (Solute carrier family 1 member 3), SLC1A3, EAAT1 GLAST GLAST1 |
| Immunogen Species   | Homo sapiens (Human)  |
| Research Area       | Others  |
| Target Names        | SLC1A3  |