




SLC22A10 Antibody

Product Code	CSB-PA723708LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63ZE4
Immunogen	Recombinant Human Solute carrier family 22 member 10 protein (37-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:100-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 10 (Organic anion transporter 5), SLC22A10, OAT5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A10

Image



Immunofluorescence staining of SH-SY5Y cells with CSB-PA723708LA01HU at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.