



SRC Antibody, Biotin conjugated

Product Code	CSB-PA022650LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12931
Immunogen	Recombinant Human Proto-oncogene tyrosine-protein kinase Src protein (2-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proto-oncogene tyrosine-protein kinase Src (EC 2.7.10.2) (Proto-oncogene c-Src) (pp60c-src) (p60-Src), SRC, SRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SRC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.