



SUMF2 Antibody, HRP conjugated

Product Code	CSB-PA839844LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBJ7
Immunogen	Recombinant Human Sulfatase-modifying factor 2 protein (26-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfatase-modifying factor 2 (C-alpha-formylglycine-generating enzyme 2), SUMF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SUMF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.