



SUMF2 Antibody, FITC conjugated

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA839844LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NBJ7 |
| Immunogen | Recombinant Human Sulfatase-modifying factor 2 protein (26-301AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Antigen Affinity purified & Affinity purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sulfatase-modifying factor 2 (C-alpha-formylglycine-generating enzyme 2), SUMF2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SUMF2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |