



TAS2R5 Antibody, Biotin conjugated

Product Code	CSB-PA885769PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYW4
Immunogen	Peptide sequence from Human Taste receptor type 2 member 5 protein (149-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 5 (T2R5), TAS2R5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAS2R5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.