




TAS2R7 Antibody

Product Code	CSB-PA868363OA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYW3
Immunogen	Peptide sequence from Human Taste receptor type 2 member 7 protein (151-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 7 (T2R7) (Taste receptor family B member 4) (TRB4), TAS2R7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAS2R7

Image



Western Blot

Positive WB detected in: Hela whole cell lysate, A549 whole cell lysate

All lanes: TAS2R7 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 37 kDa

Observed band size: 37 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.