



TBC1D3G Antibody, HRP conjugated

Product Code	CSB-PA732237LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DHY5
Immunogen	Recombinant Human TBC1 domain family member 3G protein (312-549AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 3G, TBC1D3G
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D3G