



TBC1D3G Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA732237LC11HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6DHY5 |
| Immunogen | Recombinant Human TBC1 domain family member 3G protein (312-549AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | TBC1 domain family member 3G, TBC1D3G |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | TBC1D3G |