



TMEM37 Antibody, HRP conjugated

Product Code	CSB-PA820182NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXS4
Immunogen	Peptide sequence from Human Voltage-dependent calcium channel gamma-like subunit protein (4-20AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent calcium channel gamma-like subunit (Neuronal voltage-gated calcium channel gamma-like subunit) (Transmembrane protein 37), TMEM37, PR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM37