



TMEM37 Antibody, FITC conjugated

Product Code	CSB-PA820182NC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXS4
Immunogen	Peptide sequence from Human Voltage-dependent calcium channel gamma-like subunit protein (4-20AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent calcium channel gamma-like subunit (Neuronal voltage-gated calcium channel gamma-like subunit) (Transmembrane protein 37), TMEM37, PR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM37
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.