



VPS26A Antibody, Biotin conjugated

Product Code	CSB-PA025898LD11HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75436
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 26A protein (21-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 26A (Vesicle protein sorting 26A) (hVPS26), VPS26A, VPS26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VPS26A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.