



VPS26A Antibody, Biotin conjugated

Product Code	CSB-PA025898LD11HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75436
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 26A protein (21-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 26A (Vesicle protein sorting 26A) (hVPS26), VPS26A, VPS26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VPS26A