



YAP1 Antibody, Biotin conjugated

Product Code	CSB-PA026244LD11HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46937
Immunogen	Recombinant Human Transcriptional coactivator YAP1 protein (155-504AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional coactivator YAP1 (Yes-associated protein 1) (Protein yorkie homolog) (Yes-associated protein YAP65 homolog), YAP1, YAP65
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	YAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.