



Apoa1 Antibody, Biotin conjugated

Product Code	CSB-PA001913LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04639
Immunogen	Recombinant Rat ApoA1 protein (19-259AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein A-I (Apo-AI) (ApoA-I) (Apolipoprotein A1) [Cleaved into: Proapolipoprotein A-I (ProapoA-I)], ApoA1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Apoa1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.