



CA9 Antibody, Biotin conjugated

Product Code	CSB-PA614990LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16790
Immunogen	Peptide sequence from Human Carbonic anhydrase 9 protein (38-414AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonic anhydrase 9 (EC 4.2.1.1) (Carbonate dehydratase IX) (Carbonic anhydrase IX) (CA-IX) (CAIX) (Membrane antigen MN) (P54/58N) (Renal cell carcinoma-associated antigen G250) (RCC-associated antigen G250) (pMW1), CA9, G250 MN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CA9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.