



GCG Antibody, Biotin conjugated

Product Code	CSB-PA009315LD11HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01275
Immunogen	Recombinant Human Glucagon protein (21-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucagon [Cleaved into: Glicentin, Glicentin-related polypeptide (GRPP), Oxyntomodulin (OXM) (OXY), Glucagon, Glucagon-like peptide 1 (GLP-1) (Incretin hormone), Glucagon-like peptide 1(7-37) (GLP-1(7-37)), Glucagon-like peptide 1(7-36) (GLP-1(7-36)), Glucagon-like peptide 2 (GLP-2)], GCG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GCG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.