



HLA-G Antibody, Biotin conjugated

Product Code	CSB-PA010509LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17693
Immunogen	Recombinant Human HLA class I histocompatibility antigen, alpha chain G protein (25-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class I histocompatibility antigen, alpha chain G (HLA G antigen) (MHC class I antigen G), HLA-G, HLA-6.0 HLAG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HLA-G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.