




Recombinant Human Annexin A5 (ANXA5)

Product Code	CSB-EP001846HU
Abbreviation	Recombinant Human ANXA5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08758
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AQVLRGTVTDFPGFDERADAETLRKAMKGLGTDEESILTLTSSRNAQRQEIS AAFKTLFGRDLLDDLKSELTKFEKLIVALMKPSRLYDAYELKHALKGAGTNEK VLTEIIASRTPEELRAIKQVYEEYEGSSLEDDVVGDTSGYYQRMLVLLQANRD PDAGIDEAQVEQDAQALFQAGELKWGTDEEKFITIFGTRSVSHLRKVFDDKYMTI SGFQIEETIDRETSGNLEQLLLAVVKSIRSIPAYLAETLYAMKGAGTDDHTLIRV MVSRSSEIDLFNIRKEFRKNFATSLYSMIKGDTS GDYK KALLLLCGEDD
Research Area	Cancer
Source	E.coli
Target Names	ANXA5
Protein Names	Recommended name: Annexin A5 Alternative name(s): Anchorin CII Annexin V Annexin-5 Calphobindin I Short name= CBP-I Endonexin II Lipocortin V Placental anticoagulant protein 4 Short name= PP4 Placental ant
Expression Region	2-320aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	62.8 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.