



Recombinant *Cyberlindnera jadinii* Uricase

Product Code	CSB-MP025648EPX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P78609
Product Type	Recombinant Protein
Immunogen Species	<i>Cyberlindnera jadinii</i> (Torula yeast) (<i>Pichia jadinii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSTTSSSTY GKDNVFKLV KKDPQNPQQ EVMEATVTCL LEGGFDTSTY EADNSSIVPT DTVKNTILVL AKTTEIWPIE RFAAKLATHF VEKYSHVSGV SVKIVQDRWV KYAVDGKPHD HSFIEGGEK RITDLYYKRS GDYKLSSAIK DLTVLKSTGS MFYGYNKCDF TTLQPTTDRI LSTDVDATWV WDNKKIGSVY DIAKAADKGI FDNVYNQARE ITLTTFALFN SPSVQATMFN MATQILEKAC SVYSVSYALP NKHYFLIDLK WKGLENDNEL FYSPHPNGL IKCTVVRKEK TKL
Source	Mammalian cell
Protein Names	Recommended name: Uricase EC= 1.7.3.3 Alternative name(s): Urate oxidase
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.