




Recombinant Corynephage omega Diphtheria toxin, partial

Product Code	CSB-EP360420CQS
Abbreviation	Recombinant Corynephage omega Diphtheria toxin protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00587
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Corynephage omega
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	GADDVVDSSKSFVMENFSSYHGTKPGYVDSIQKGIQKPKSGTQGNYYDDDWK GFYSTDNKYDAAGYSVDNENPLSGKAGGVVKVYTPGLTKVLALKVDNAETIKK ELGLSLTEPLMEQVGTEEFIKRFGDGASRVVLSLPFAEGSSSVEYINNWEQAK ALSVELEINFETRGRGQDAMYEYMAQACAGNRVRR
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Diphtheria toxin Short name= DTA Alternative name(s): NAD(+)-diphthamide ADP-ribosyltransferase EC= 2.4.2.36 Cleaved into the following 2 chains: 1. Diphtheria toxin fragment A 2. Diphtheria toxin fragment B
Expression Region	26-218aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	28.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.