




Recombinant Human Flap endonuclease 1 (FEN1)

Product Code	CSB-EP008585HU
Abbreviation	Recombinant Human FEN1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39748
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MGIQGLAKLIADVAPSAIRENDIKSYFGRKVAIDASMSIYQFLIAVRQGGDVLQN EEGETTSHLMGMFYRTIRMMENGIKPVYVFDGKPPQLKSGELAKRSERRAEA EKQLQQAQAAGAEQEVEKFTKRLVKVTKQHNDCKHLLSLMGIPYLDAPSEAE ASCAALVKAGKVYAAATEDMDCLTFGSPVLMRHLTASEAKKLP IQEFHLSRILQ ELGLNQE QFVDLCILLGSDYCESIRGIGPKRAVDLIQKHKSIEEIVRRLDPNKYP VPENWLHKEAHQLFLEPEVLDPEVELKWSEPNEEELIKFMCGEKQFSEERIR SGVKRLSKSRQGSTQGRLLDDFFKVTGSLSSAKRKEPEPKGSTKKKAKTGAAG KFKRGK
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): DNase IV Flap structure-specific endonuclease 1 Maturation factor 1 Short name= MF1 Short name= hFEN-1
Expression Region	1-380aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	50.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.