




Recombinant Human Fc receptor-like protein 5 (FCRL5), partial

Product Code	CSB-MP853494HU
Abbreviation	Recombinant Human FCRL5 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96RD9
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QFARTPRPIIFLQPPWTTVFQGERVTLTCKGFRFYSPQKTKWYHRYLGKEILRE TPDNILEVQESGEYRCQAQGSPLSSPVHLDLFSSASLILQAPLSVFEGDSVVLRC RAKAEVTLNNTIYKNDNVLAFLNKRTDFHIPHACKDNGAYRCTGYKESCCPV SSNTVKIQVQEPFTRPVLRASSFPISGNPVTLCETQLSLERSDVPLRFRFFR DDQTLGLGWSLSPNFQITAMWSKDSGFYWCKAATMPYSVISDSPRSWIQVQI PASHPVLTLSPKALNFEGTKVTLHCETQEDSLRTLYRFYHEGVPLRHKSVCRC ERGASISFSLTTENSGNYYCTADNGLGAKPSKAVSLSVTPVSHPVNLSSPE DLIFEGAKVTLHCEAQRGSLPILYQFHHEGAALERRSANSAGGVAISFSLTAEH SGNYYCTADNGFGPQRSKAVSLSVTPVSHPVLTLSAEALTFEGATVTLHCE VQRGSPQILYQFYHEDMPLWSSSTPSVGRVSFSFSLTEGHSGNYYCTADNGF GPQRSEVVSFLFVTPVSRPILTLRVPRQAQAVVGDILLEHCEAPRGSPPILYWYF HEDVTLGSSSAPSGGEASFNLSTAEHSGNYSCEANGLVAQHSdTISLSVIV PVSRPILTFRAPRAQAVVGDILLEHCEALRGSSPILYWYFHEVTLGKISAPSG GGASFNLSTTEHSGIYSCEADNGLEAQRSEMVTLKVAVPVSRPVLTLRAPGT HAAVGDILLEHCEALRGSPILYRFFHEDVTLGNRSSPSGGASLNLSTAEHS GNYSCEADNGLGAQRSETVTLYITGLTANR
Research Area	Immunology
Source	Mammalian cell
Target Names	FCRL5
Protein Names	Recommended name: Fc receptor-like protein 5 Short name= FcR-like protein 5 Short name= FcRL5 Alternative name(s): BXMAS1 Fc receptor homolog 5 Short name= FcRH5 Immune receptor translocation-associated protein 2 C
Expression Region	16-844aa



Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	92.9 kDa
Protein Length	Partial

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-----------------------	---

Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
-------------------	---