



Goat Anti-Human IgG, Fcγ fragment specific; HRP conjugated

Product Code	CSB-PA160506
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: ELISA:1:5000-1:100000, WB:1:5000-1:50000
Form	Liquid
Conjugate	HRP
Storage Buffer	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Isotype	IgG
Clonality	Polyclonal
Alias	IgG, Fcγ fragment specific
Product Type	Secondary Antibody
Immunogen Species	Human