





# Recombinant Severe acute respiratory syndrome coronavirus 2 Replicase polyprotein 1ab (rep), partial

<b>Product Code</b>	CSB-MP3388GND
<b>Abbreviation</b>	Recombinant SARS-CoV-2 Replicase polyprotein 1ab protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0DTD1
<b>Form</b>	Lyophilized powder
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE. ≥ 95% as determined by SEC-HPLC.
<b>Sequence</b>	NNELSPVALRQMSCAAGTTQACTDDNALAYYNTTKGGRFVLALLSDLQLK WARFPKSDGTGTIYTELEPPCRFVTDTPKGPKVKYLYFIKGLNNLNRMVLGS LAATVRLQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rep
<b>Protein Names</b>	Non-structural protein 9
<b>Expression Region</b>	4141-4253aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal hFc1-Flag-tagged
<b>Mol. Weight</b>	44.5 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of rep was greater than 95% as determined by SEC-HPLC

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.