



# Goat Anti-Human IgG, Fc $\gamma$ fragment specific Antibody

<b>Product Code</b>	CSB-PA793692
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Goat
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Secondary Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.