



Goat Anti-Human IgG, Fcγ fragment specific Antibody

Product Code	CSB-PA793692
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Isotype	IgG
Clonality	Polyclonal
Product Type	Secondary Antibody
Immunogen Species	Homo sapiens (Human)