



# N Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA33254C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0DTC9
<b>Immunogen</b>	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa) (95%, Protein G purified)
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Research Area</b>	Microbiology